



CXXC1 Rabbit pAb

Catalog No	YP-Ab-19246
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	CXXC1 CFP1 CGBP PCCX1 PHF18
Protein Name	CpG-binding protein (CXXC-type zinc finger protein 1) (PHD finger and CXXC domain-containing protein 1)
Immunogen	Synthesized peptide derived from human CXXC1
Specificity	This antibody detects endogenous levels of CXXC1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	72kD
Cell Pathway	Nucleus speckle . Nucleus . Associated with euchromatin. During mitosis, excluded from condensed chromosomes.
Tissue Specificity	Ubiquitous.
Function	Transcriptional activator that exhibits a unique DNA binding specificity for CpG unmethylated motifs with a preference for CpGG.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images