



SNRPG Rabbit pAb

Catalog No	YP-Ab-19251
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	SNRPG PBSCG
Protein Name	Small nuclear ribonucleoprotein G (snRNP-G) (Sm protein G) (Sm-G) (SmG)
Immunogen	Synthesized peptide derived from human SNRPG
Specificity	This antibody detects endogenous levels of SNRPG at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	8kD
Cell Pathway	Cytoplasm, cytosol . Nucleus . SMN-mediated assembly into core snRNPs occurs in the cytosol before SMN-mediated transport to the nucleus to be included in spliceosomes. .
Tissue Specificity	
Function	Plays a role in pre-mRNA splicing as a core component of the spliceosomal U1, U2, U4 and U5 small nuclear ribonucleoproteins (snRNPs), the building blocks of the spliceosome . Component of both the pre-catalytic spliceosome B complex and activated spliceosome C complexes . Is also a component of the minor U12 spliceosome . As part of the U7 snRNP it is involved in histone 3'-end processing .
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images