



METTL2A Mouse mAb

Catalog No	YP-mAb-19250
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	
Protein Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	METTL2; METTL2A
Observed Band	
Calculated Molecular Weight	44kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Enables tRNA (cytosine-3-)-methyltransferase activity. Involved in tRNA methylation.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images