



RS26 mouse mAb

Catalog No	YP-mAb-09185
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPS26
Protein Name	40S ribosomal protein S26
Immunogen	Synthesized peptide derived from human RS26 AA range: 59-109
Specificity	This antibody detects endogenous levels of human RS26
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	13kD
Synonyms	
Cell Pathway	Cytoplasm, cytosol . Cytoplasm . Rough endoplasmic reticulum . Detected on cytosolic polysomes (PubMed:25957688). Detected in ribosomes that are associated with the rough endoplasmic reticulum (By similarity). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

