



DNAJC1 Rabbit pAb

Catalog No	YP-Ab-19275
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	DNAJC1 HTJ1
Protein Name	DnaJ homolog subfamily C member 1 (DnaJ protein homolog MTJ1)
Immunogen	Synthesized peptide derived from human DNAJC1
Specificity	This antibody detects endogenous levels of DNAJC1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	61kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type I membrane protein . Nucleus membrane; Single-pass type I membrane protein. Microsome membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	May modulate protein synthesis.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images