



PCIF1 Rabbit pAb

Catalog No	YP-Ab-19291
Isotype	IgG
Reactivity	Human,Mouse
Applications	
Gene Name	PCIF1 C20orf67
Protein Name	Phosphorylated CTD-interacting factor 1
Immunogen	Synthesized peptide derived from human PCIF1
Specificity	This antibody detects endogenous levels of PCIF1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	77kD
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitous.
Function	Cap-specific adenosine methyltransferase that catalyzes formation of N(6),2'-O-dimethyladenosine cap (m6A(m)) by methylating the adenosine at the second transcribed position of capped mRNAs . Recruited to the early elongation complex of RNA polymerase II (RNAPII) via interaction with POLR2A and mediates formation of m6A(m) co-transcriptionally .
Background	
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images