



FCER1A mouse mAb

Catalog No	YP-mAb-12578
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	FCER1A FCE1A
Protein Name	FCER1A
Immunogen	Synthesized peptide derived from human FCER1A
Specificity	This antibody detects endogenous levels of FCER1A at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.70% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	High affinity immunoglobulin epsilon receptor subunit alpha (Fc-epsilon RI-alpha) (FcERI) (IgE Fc receptor subunit alpha)
Calculated Molecular Weight	28kD
Cell Pathway	Cell membrane; Single-pass type I membrane protein.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

