



HELQ rabbit pAb

Catalog No	YP-Ab-09459
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	HELQ HEL308
Protein Name	Helicase POLQ-like (EC 3.6.4.12) (Mus308-like helicase) (POLQ-like helicase)
Immunogen	Synthesized peptide derived from human HELQ AA range: 830-880
Specificity	This antibody detects endogenous levels of human HELQ
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Calculated Molecular Weight	121kD
Cell Pathway	Nucleus . Chromosome . Localizes to sites of DNA damage; localizes to damaged replication forks. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

