



NBEAL2 Rabbit mAb

Catalog No	YP-rAb-17696
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	NBEAL2;KIAA0540;UNQ253/PRO290
Protein Name	Neurobeachin-like protein 2
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-5000
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15° C to -25° C/1 year(Do not lower than -25° C)
Synonyms	NBEAL2 ; KIAA0540 ; UNQ253/PRO290 ; Neurobeachin-like protein 2 ;
Observed Band	303kD
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum.
Tissue Specificity	Expressed in megakaryocytes. {ECO:0000269 PubMed:21765411}.
Function	Probably involved in thrombopoiesis. Plays a role in the development or secretion of alpha-granules, that contain several growth factors important for platelet biogenesis. {ECO:0000269 PubMed:21765411, ECO:0000269 PubMed:21765412}.
Background	The protein encoded by this gene contains a beige and Chediak-Higashi (BEACH) domain and multiple WD40 domains, and may play a role in megakaryocyte alpha-granule biogenesis. Mutations in this gene are a cause of gray platelet syndrome. [provided by RefSeq, Dec 2011]
matters needing attention	Avoid repeated freezing and thawing!





Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

HepG2

kDa
200-
140-
100-
70-
55-
40-
35-
25-
20-
15-
10-



NBEAL2

Western Blot analysis of HepG2 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-NBEAL2 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.





杭州臻优品生物科技有限公司

热销产品:

蛋白、一抗、抗体对、ELISA试剂盒、生化试剂盒
CCK8试剂盒、QPCR检测试剂盒

检测服务:

ELISA检测及定制服务 | 生化检测 | PCR、QPCR检测 | WB检测
ICO-IP检测 | 切片 | 染色 | 免疫组化 | 免疫荧光 | 透射电镜全套
| 宏基因组、转录组、基因组、蛋白组、代谢组测序



关注官网



关注客服