



# LAMP1 Rabbit mAb (AbFluor 488)

<b>Catalog No</b>	YP-rAb-18126
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human
<b>Applications</b>	IF
<b>Gene Name</b>	LAMP1
<b>Protein Name</b>	Lysosome-associated membrane glycoprotein 1
<b>Purification Process</b>	Protein A
<b>Specificity</b>	Endogenous
<b>Formulation</b>	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
<b>Source</b>	Monoclonal, Rabbit,IgG
<b>Dilution</b>	IF 1:200-1000
<b>Concentration</b>	0.5 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-15° C to -25° C/1 year(Do not lower than -25° C)
<b>Synonyms</b>	LAMP1 ; Lysosome-associated membrane glycoprotein 1 ; LAMP-1 ; Lysosome-associated membrane protein 1 ; CD107 antigen-like family member A ; CD107a
<b>Observed Band</b>	100kD
<b>Calculated Molecular Weight</b>	45kD
<b>Cell Pathway</b>	Membrane
<b>Tissue Specificity</b>	Brain,Liver,Pancreas,Placenta,Plasma,Testis,
<b>Function</b>	Presents carbohydrate ligands to selectins. Also implicated in tumor cell metastasis.,PTM:O- and N-glycosylated; some of the 18 N-linked glycans are polylactosaminoglycans.,similarity:Belongs to the LAMP family.,subcellular location:This protein shuttles between lysosomes, endosomes, and the plasma membrane.,
<b>Background</b>	The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may also play a role in tumor cell metastasis. [provided by RefSeq, Jul 2008],





**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

杭州臻优品生物科技有限公司

热销产品:

蛋白、一抗、抗体对、ELISA试剂盒、生化试剂盒  
CCK8试剂盒、QPCR检测试剂盒

检测服务:

ELISA检测及定制服务 | 生化检测 | PCR、QPCR检测 | WB检测  
ICO-IP检测 | 切片 | 染色 | 免疫组化 | 免疫荧光 | 透射电镜全套  
| 宏基因组、转录组、基因组、蛋白组、代谢组测序



关注官网



关注客服