



# Melan-A (ABT94R) Rabbit mAb (Ready to Use)

<b>Catalog No</b>	YP-rAb-18142
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human
<b>Applications</b>	IHC
<b>Gene Name</b>	MLANA
<b>Protein Name</b>	Mart-1/Melan-A
<b>Purification Process</b>	Protein A
<b>Specificity</b>	This antibody detects endogenous levels of Melan-A
<b>Formulation</b>	The prediluted ready-to-use antibody is diluted in phosphate buffer saline containing stabilizing protein and 0.05% Proclin 300
<b>Source</b>	Monoclonal, Rabbit,IgG
<b>Dilution</b>	Ready to use for IHC Note: For IHC, we suggest antigen retrieval with TE buffer pH 9.0
<b>Concentration</b>	0.5 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	2° C to 8° C/1 year,Ship by ice bag
<b>Synonyms</b>	Antigen LB39 AA ; Antigen LB39-AA ; Antigen SK29 AA ; Antigen SK29-AA ; MAR1_HUMAN ; MART 1 ; MART-1 ; MART1 ; Melan A ; Melan A protein ; Melanoma antigen recognized by T cells 1 ; Melanoma antigen recognized by T-cells 1 ; MLAN A ; MLANA ; OTTHUMP00000021036 ; OTTHUMP00000021037 ; OTTHUMP00000021038 ; Protein Melan-A
<b>Observed Band</b>	
<b>Calculated Molecular Weight</b>	
<b>Cell Pathway</b>	Cytoplasmic
<b>Tissue Specificity</b>	Expression is restricted to melanoma and melanocyte cell lines and retina.
<b>Function</b>	tissue specificity:Expression is restricted to melanoma and melanocyte cell lines and retina.,
<b>Background</b>	Tissue specificity:Expression is restricted to melanoma and melanocyte cell lines and retina.,





**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**杭州臻优品生物科技有限公司**

**热销产品:**

蛋白、一抗、抗体对、ELISA试剂盒、生化试剂盒  
CCK8试剂盒、QPCR检测试剂盒

**检测服务:**

ELISA检测及定制服务 | 生化检测 | PCR、QPCR检测 | WB检测  
ICO-IP检测 | 切片 | 染色 | 免疫组化 | 免疫荧光 | 透射电镜全套  
| 宏基因组、转录组、基因组、蛋白组、代谢组测序



关注官网



关注客服