



# CHCH3 rabbit pAb

<b>Catalog No</b>	YP-Ab-09824
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	CHCHD3 MINOS3
<b>Protein Name</b>	Coiled-coil-helix-coiled-coil-helix domain-containing protein 3, mitochondrial
<b>Immunogen</b>	Synthesized peptide derived from human CHCH3 AA range: 144-194
<b>Specificity</b>	This antibody detects endogenous levels of human CHCH3
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:1000-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	25kD
<b>Cell Pathway</b>	Mitochondrion inner membrane ; Lipid-anchor ; Intermembrane side . Cytoplasm . Nucleus . Mitochondrion .
<b>Tissue Specificity</b>	Detected at low levels in brain, placenta, lung, liver, kidney and pancreas with increased levels in heart and skeletal muscle. Higher expression in primary lung cancers than in normal lung tissue.
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

