



DAND5 Rabbit pAb

Catalog No	YP-Ab-19310
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	DAND5 CER2 CKTSF1B3 GREM3 SP1
Protein Name	DAN domain family member 5 (Cerberus-like protein 2) (Cerl-2) (Cysteine knot superfamily 1, BMP antagonist 3) (Gremlin-3)
Immunogen	Synthesized peptide derived from human DAND5
Specificity	This antibody detects endogenous levels of DAND5 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	21kD
Cell Pathway	Secreted .
Tissue Specificity	
Function	Seems to play a role in the correct specification of the left-right axis. May antagonize NODAL and BMP4 signaling. Cystine knot-containing proteins play important roles during development, organogenesis, tissue growth and differentiation (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images