



# PUF60 Rabbit pAb

|                                    |   |
|------------------------------------|---|
| <b>Catalog No</b>                  | YP-Ab-19317   |
| <b>Isotype</b>                     | IgG   |
| <b>Reactivity</b>                  | Human;Mouse;Rat   |
| <b>Applications</b>                | WB  |
| <b>Gene Name</b>                   | PUF60 FIR ROBPI SIAHBP1   |
| <b>Protein Name</b>                | Poly(U)-binding-splicing factor PUF60 (60 kDa poly(U)-binding-splicing factor) (FUSE-binding protein-interacting repressor) (FBP-interacting repressor) (Ro-binding protein 1) (RoBP1) (Siah-binding protein 1) (Siah-BP1)  |
| <b>Immunogen</b>                   | Synthesized peptide derived from human PUF60  |
| <b>Specificity</b>                 | This antibody detects endogenous levels of PUF60 at Human, Mouse,Rat  |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| <b>Source</b>                      | Polyclonal, Rabbit,IgG  |
| <b>Purification</b>                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.   |
| <b>Dilution</b>                    | WB 1:500-2000   |
| <b>Concentration</b>               | 1 mg/ml   |
| <b>Purity</b>                      | ≥90%  |
| <b>Storage Stability</b>           | -20°C/1 year  |
| <b>Synonyms</b>                    |   |
| <b>Observed Band</b>               |   |
| <b>Calculated Molecular Weight</b> | 61kD  |
| <b>Cell Pathway</b>                | Nucleus . Colocalizes partially with RO60. .  |
| <b>Tissue Specificity</b>          | Isoform 2 is expressed in colonic epithelium and colorectal epithelium cancer (at protein level). Isoform 6 is expressed in colorectal epithelial cancer but below detection level in colonic epithelium. Expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney, pancreas, spleen, thymus, prostate, testis, ovary, small intestine, colon and peripheral blood leukocytes.   |
| <b>Function</b>                    | DNA- and RNA-binding protein, involved in several nuclear processes such as pre-mRNA splicing, apoptosis and transcription regulation. In association with FUBP1 regulates MYC transcription at the P2 promoter through the core-TFIID basal transcription factor. Acts as a transcriptional repressor through the core-TFIID basal transcription factor. Represses FUBP1-induced transcriptional activation but not basal transcription. Decreases ERCC3 helicase activity. Does not repress TFIID-mediated transcription in xeroderma pigmentosum complementation group B (XPB) cells. Is also involved in pre-mRNA splicing. Promotes splicing of an intron with weak 3'-splice site and pyrimidine tract in a |



cooperative manner with U2AF2. Involved in apoptosis induction when overexpressed in HeLa cells. Isoform 6 failed to repress MYC transcription and inhibited FIR-induced apoptosis in colorectal cancer. Isoform 6 may contribute to tumor progression by enabling increased MYC expression and greater resistance to apoptosis in tumors than in normal cells. Modulates alternative splicing of several mRNAs. Binds to relaxed DNA of active promoter regions. Binds to the pyrimidine tract and 3'-splice site regions of pre-mRNA; binding is enhanced in presence of U2AF2. Binds to Y5 RNA in association with RO60. Binds to poly(U) RNA.

## Background

### matters needing attention

Avoid repeated freezing and thawing!

### Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images