



CPPED mouse mAb

Catalog No	YP-mAb-12144
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	CPPED1 CSTP1
Protein Name	CPPED
Immunogen	Synthesized peptide derived from human CPPED AA range: 68-118
Specificity	This antibody detects endogenous levels of CPPED at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	35kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed in subcutaneous adipose tissue.
Function	cofactor: Binds 2 divalent metal cations.,induction: Transactivated by the large envelope protein of the hepatitis B virus (HBV).,similarity: Belongs to the metallophosphoesterase superfamily. CPPED1 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using CPPED mouse mAb

