



MT1M Mouse mAb

Catalog No	YP-mAb-18423
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	MT1M MT1K
Protein Name	Metallothionein-1M (MT-1M) (Metallothionein-IM) (MT-IM)
Immunogen	Synthesized peptide derived from human MT1M
Specificity	This antibody detects endogenous levels of MT1M at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	7kD
Cell Pathway	
Tissue Specificity	
Function	Metallothioneins have a high content of cysteine residues that bind various heavy metals; these proteins are transcriptionally regulated by both heavy metals and glucocorticoids.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

