



Actin beta rabbit pAb(Yeast)

Catalog No	YP-Ab-10305
Isotype	IgG
Reactivity	Yeast
Applications	WB
Gene Name	
Protein Name	Actin beta
Immunogen	Purified recombinant Full length Actin beta (Caenorhabditis elegans)
Specificity	This antibody detects endogenous levels of Actin beta (Caenorhabditis elegans)
Formulation	Liquid in PBS containing 50% glycerol and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Beta Actin,βActin
Observed Band	42kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of Yeast lysis, using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000

