




ERLN1 mouse mAb

Catalog No	YP-mAb-12104
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ERLN1 C10orf69 KE04 KEO4 SPFH1
Protein Name	ERLN1
Immunogen	Synthesized peptide derived from human ERLN1 AA range: 185-235
Specificity	This antibody detects endogenous levels of ERLN1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	38kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type II membrane protein . Associated with lipid raft-like domains of the endoplasmic reticulum membrane.
Tissue Specificity	Expressed in heart, placenta, liver, kidney, pancreas, prostate, testis, ovary and small intestine.
Function	similarity:Belongs to the band 7/mec-2 family.,subcellular location:Associated with lipid raft-like domains of the endoplasmic reticulum membrane.,tissue specificity:Expressed in heart, placenta, liver, kidney, pancreas, prostate, testis, ovary and small intestin.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ERLN1 mouse mAb