



SHANK1 Rabbit pAb

Catalog No	YP-Ab-19339
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SHANK1
Protein Name	SH3 and multiple ankyrin repeat domains protein 1 (Shank1) (Somatostatin receptor-interacting protein) (SSTR-interacting protein) (SSTRIP)
Immunogen	Synthesized peptide derived from human SHANK1. AA range:372-472
Specificity	This antibody detects endogenous levels of SHANK1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	238kD
Cell Pathway	Cytoplasm . Cell junction, synapse, postsynaptic density . Cell junction, synapse . Colocalizes with alpha-latrotoxin receptor 1. .
Tissue Specificity	Expressed in brain particularly in the amygdala, hippocampus, substantia nigra and thalamus. Isoform 2 seems to be expressed ubiquitously.
Function	Seems to be an adapter protein in the postsynaptic density (PSD) of excitatory synapses that interconnects receptors of the postsynaptic membrane including NMDA-type and metabotropic glutamate receptors via complexes with GKAP/PSD-95 and Homer, respectively, and the actin-based cytoskeleton. Plays a role in the structural and functional organization of the dendritic spine and synaptic junction.
Background	
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images