



PGP Rabbit pAb

Catalog No	YP-Ab-19355
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	PGP
Protein Name	Phosphoglycolate phosphatase (PGP) (PGPase) (EC 3.1.3.18)
Immunogen	Synthesized peptide derived from human PGP. AA range:1-99
Specificity	This antibody detects endogenous levels of PGP at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	35kD
Cell Pathway	
Tissue Specificity	Detected in all tissues including red cells, lymphocytes and cultured fibroblasts (at protein level). The highest activities occur in skeletal muscle and cardiac muscle.
Function	Glycerol-3-phosphate phosphatase hydrolyzing glycerol-3-phosphate into glycerol. Thereby, regulates the cellular levels of glycerol-3-phosphate a metabolic intermediate of glucose, lipid and energy metabolism. Was also shown to have a 2-phosphoglycolate phosphatase activity and a tyrosine-protein phosphatase activity. However, their physiological relevance is unclear. In vitro, has also a phosphatase activity toward ADP, ATP, GDP and GTP (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images