




ISPD mouse mAb

Catalog No	YP-mAb-12004
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ISPD
Protein Name	ISPD
Immunogen	Synthesized peptide derived from human ISPD AA range: 301-351
Specificity	This antibody detects endogenous levels of ISPD at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	50kD
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Ubiquitously expressed, with high expression in brain.
Function	similarity:Belongs to the ispD family.,
Background	This gene encodes a 2-C-methyl-D-erythritol 4-phosphate cytidyltransferase-like protein. Mutations in this gene are the cause of Walker-Warburg syndrome. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2012],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ISPD mouse mAb