



LYOX rabbit pAb

Catalog No	YP-Ab-09451
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LOX
Protein Name	Protein-lysine 6-oxidase (EC 1.4.3.13) (Lysyl oxidase)
Immunogen	Synthesized peptide derived from human LYOX
Specificity	This antibody detects endogenous levels of human LYOX
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	46kD
Cell Pathway	Secreted . Secreted, extracellular space.
Tissue Specificity	Heart, placenta, skeletal muscle, kidney, lung and pancreas.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

