



DI3L2 mouse mAb

Catalog No	YP-mAb-09907
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	DIS3L2 FAM6A
Protein Name	DIS3-like exonuclease 2 (EC 3.1.13.-)
Immunogen	Synthesized peptide derived from human DI3L2 AA range: 306-356
Specificity	This antibody detects endogenous levels of human DI3L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	97kD
Cell Pathway	Cytoplasm . Cytoplasm, P-body.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images