



ZFHX4 mouse mAb

Catalog No	YP-mAb-12566
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ZFHX4
Protein Name	ZFHX4
Immunogen	Synthesized peptide derived from human ZFHX4
Specificity	This antibody detects endogenous levels of ZFHX4 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.57% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zinc finger homeobox protein 4 (Zinc finger homeodomain protein 4) (ZFH-4)
Observed Band	392kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in brain, skeletal muscle and liver. Very low expression in stomach.
Function	disease:A chromosomal aberration involving [ZFHX4] is found in one patient with ptosis. Translocation t(1;8)(p34.3;q21.12).,function:May play a role in neural and muscle differentiation (By similarity). May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 20 C2H2-type zinc fingers.,similarity:Contains 4 homeobox DNA-binding domains.,tissue specificity:Expressed in brain, skeletal muscle and liver. Very low expression in stomach.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images