



RPL21 Mouse mAb

Catalog No	YP-mAb-19369
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	RPL21
Protein Name	60S ribosomal protein L21
Immunogen	Synthesized peptide derived from human RPL21. AA range:1-130
Specificity	This antibody detects endogenous levels of RPL21 at Human,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	18kD
Cell Pathway	Cytoplasm, cytosol . Cytoplasm . Endoplasmic reticulum . Detected on cytosolic polysomes (PubMed:25957688). Detected in ribosomes that are associated with the rough endoplasmic reticulum (By similarity). .
Tissue Specificity	
Function	Component of the large ribosomal subunit.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images