



SRGAP3 Rabbit pAb

Catalog No	YP-Ab-19383
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	SRGAP3 ARHGAP14 KIAA0411 KIAA1156 MEGAP SRGAP2
Protein Name	SLIT-ROBO Rho GTPase-activating protein 3 (srGAP3) (Mental disorder-associated GAP) (Rho GTPase-activating protein 14) (WAVE-associated Rac GTPase-activating protein) (WRP)
Immunogen	Synthesized peptide derived from human SRGAP3. AA range:59-129
Specificity	This antibody detects endogenous levels of SRGAP3 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	121kD
Cell Pathway	
Tissue Specificity	Highly expressed in adult and fetal brain. Expressed at low levels in kidney. Isoform 3 is expressed in the kidney but is absent in the brain.
Function	GTPase-activating protein for RAC1 and perhaps Cdc42 , but not for RhoA small GTPase. May attenuate RAC1 signaling in neurons.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images