



DCAF8 rabbit pAb

| | |
|----------------------------------|---|
| Catalog No | YP-Ab-09593 |
| Isotype | IgG |
| Reactivity | Human;Mouse;Rat |
| Applications | WB |
| Gene Name | DCAF8 H326 WDR42A |
| Protein Name | DDB1- and CUL4-associated factor 8 (WD repeat-containing protein 42A) |
| Immunogen | Synthesized peptide derived from human DCAF8 AA range: 65-115 |
| Specificity | This antibody detects endogenous levels of human DCAF8 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Source | Polyclonal, Rabbit,IgG |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Dilution | WB 1:1000-2000 |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | |
| Observed Band | 66kD |
| Cell Pathway | Nucleus . Cytoplasm . It shuttles between the nucleus and the cytoplasm. Nuclear import is mediated by KPNA1 and KPNB1 under the regulation of nuclear GTPase RAN. Nuclear export to the cytoplasm is XPO1 dependent. . |
| Tissue Specificity | |
| Function | |
| Background | |
| matters needing attention | Avoid repeated freezing and thawing! |
| Usage suggestions | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



Products Images