



## ZBTB24 mouse mAb

<b>Catalog No</b>	YP-mAb-19395
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human,Mouse,Rat
<b>Applications</b>	WB
<b>Gene Name</b>	ZBTB24 KIAA0441 ZNF450
<b>Protein Name</b>	Zinc finger and BTB domain-containing protein 24 (Zinc finger protein 450)
<b>Immunogen</b>	Synthesized peptide derived from human ZBTB24. AA range:429-499
<b>Specificity</b>	This antibody detects endogenous levels of ZBTB24 at Human, Mouse,Rat
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Monoclonal,Mouse,IgG
<b>Purification</b>	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Calculated Molecular Weight</b>	77kD
<b>Cell Pathway</b>	Nucleus .
<b>Tissue Specificity</b>	Widely expressed , with highest levels in naive B-cells.
<b>Function</b>	May be involved in BMP2-induced transcription.
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images