




PRR16 mouse mAb

Catalog No	YP-mAb-08713
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PRR16
Protein Name	PRR16
Immunogen	Synthesized peptide derived from human PRR16 AA range: 65-115
Specificity	This antibody detects endogenous levels of PRR16 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	33kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PRR16 mouse mAb