



TIMM17A Rabbit pAb

Catalog No	YP-Ab-19405
Isotype	IgG
Reactivity	Human, Mouse, Rat
Applications	WB
Gene Name	
Protein Name	
Immunogen	Recombinant protein (or fragment). This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Polyclonal, Rabbit, IgG
Purification	Affinity purification
Dilution	WB 1:500 - 1:2000 IF/ICC 1:50 - 1:200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TIM17; TIM17A; TIMM17A
Observed Band	18kDa
Calculated Molecular Weight	18kDa
Cell Pathway	Mitochondrion inner membrane, Multi-pass membrane protein.
Tissue Specificity	
Function	
Background	Predicted to contribute to protein transmembrane transporter activity. Predicted to be involved in protein import into mitochondrial matrix. Located in mitochondrial inner membrane and nucleoplasm. Is integral component of mitochondrial inner membrane. Part of TIM23 mitochondrial import inner membrane translocase complex.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images