




METTL14 (PT1750R)Rabbit mAb

Catalog No	YP-rAb-18535
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB、 IF、 IP、 ELISA
Gene Name	METTL14 KIAA1627
Protein Name	MET14
Immunogen	The specific immunogen used to produce this antibody is proprietary information.
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Purification	Recombinant Antibody expressed in animal component-free (ACF) media, purified via Protein A affinity chromatography.
Dilution	WB 1:1000-1:5000; IF 1:200-1:1000; ELISA 1:5000-1:20000; IP 1:50-1:200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	60kD
Calculated Molecular Weight	52kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	Function:Probable methyltransferase. ,similarity:Belongs to the MT-A70-like family. ,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-METTL14 (PT1750R) antibody. The HRP-conjugated Goat anti-Rabbit IgG (H + L) antibody was used to detect the antibody.

Lane 1: Raji
Lane 2: NIH-3T3
Lane 3: Rat testis

Predicted band size: 52kDa
Observed band size: 60kDa