



PALD mouse mAb

Catalog No	YP-mAb-09744
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PALD1 KIAA1274 PALD
Protein Name	Paladin
Immunogen	Synthesized peptide derived from human PALD AA range: 266-316
Specificity	This antibody detects endogenous levels of human PALD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	94kD
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Expressed in endothelial cells, and in certain larger vessels, in mural cells. In the brain, possibly expressed in microglia. Expressed in peripheral blood mononuclear cells (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images