



UBL3 mouse mAb

Catalog No	YP-mAb-19417
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	UBL3 PNSC1
Protein Name	Ubiquitin-like protein 3 (Membrane-anchored ubiquitin-fold protein) (HsMUB) (MUB) (Protein HCG-1)
Immunogen	Synthesized peptide derived from human UBL3. AA range:1-130
Specificity	This antibody detects endogenous levels of UBL3 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	13kD
Cell Pathway	Cell membrane ; Lipid-anchor .
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images