



DYRK2 Monoclonal Antibody

Catalog No	YP-mAb-06733
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DYRK2
Protein Name	Dual specificity tyrosine-phosphorylation-regulated kinase 2 (EC 2.7.12.1)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	DYRK2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm. Nucleus. Translocates into the nucleus following DNA damage.
Tissue Specificity	Testis, after the onset of spermatogenesis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images