



VARs2 Rabbit pAb

Catalog No	YP-Ab-19427
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	VARs2 KIAA1885 VARs2L VARSL
Immunogen	Synthesized peptide derived from human VARs2. AA range:79-179
Specificity	This antibody detects endogenous levels of VARs2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VARs2 KIAA1885 VARs2L VARSL;Valine--tRNA ligase, mitochondrial (Valyl-tRNA synthetase) (ValRS) (Valyl-tRNA synthetase-like)
Observed Band	
Calculated Molecular Weight	117kD
Cell Pathway	Mitochondrion .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

