



RIMBP2 mouse mAb

Catalog No	YP-mAb-19432
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	RIMBP2 KIAA0318 RBP2
Immunogen	Synthesized peptide derived from human RIMBP2. AA range:377-477
Specificity	This antibody detects endogenous levels of RIMBP2 at Human,Rat, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	RIMBP2 KIAA0318 RBP2;RIMS-binding protein 2 (RIM-BP2)
Observed Band	
Calculated Molecular Weight	116kD
Cell Pathway	Cell membrane . Cell junction , synapse . Synaptic plasma membrane. .
Tissue Specificity	
Function	Plays a role in the synaptic transmission as bifunctional linker that interacts simultaneously with RIMS1 , RIMS2 , CACNA1D and CACNA1B.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

