



ELAV3 rabbit pAb

Catalog No	YP-Ab-10048
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ELAVL3 HUC PLE21
Protein Name	ELAV-like protein 3 (Hu-antigen C) (HuC) (Paraneoplastic cerebellar degeneration-associated antigen) (Paraneoplastic limbic encephalitis antigen 21)
Immunogen	Synthesized peptide derived from human ELAV3 AA range: 110-160
Specificity	This antibody detects endogenous levels of human ELAV3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	40kD
Cell Pathway	
Tissue Specificity	Brain specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

