



HSC70 Rabbit mAb

Catalog No	YP-rAb-17003
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB,IHC,IF,IP,ELISA
Gene Name	HSPA8 HSC70 HSP73 HSPA10
Protein Name	Heat shock cognate 71 kDa protein
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	IHC 1:200-1:1000; WB 1:2000-1:10000; IF 1:200-1:1000; ELISA 1:5000-1:20000; IP 1:50-1:200; Note: For IHC, we suggest antigen retrieval with TE buffer pH 9.0
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	HSC70;Heat shock cognate 71 kDa protein;HSPA8 ; HSC70 ; HSP73 ; HSPA10 ; Heat shock cognate 71 kDa protein ; Heat shock 70 kDa protein 8;HSPA8 HSC70 HSP73 HSPA10
Observed Band	71kD
Calculated Molecular Weight	71kD
Cell Pathway	Cytoplasm. Melanosome. Nucleus, nucleolus. Cell membrane. Localized in cytoplasmic mRNP granules containing untranslated mRNAs. Translocates rapidly from the cytoplasm to the nuclei, and especially to the nucleoli, upon heat shock.
Tissue Specificity	Ubiquitous.
Function	Chaperone. Isoform 2 may function as an endogenous inhibitory regulator of HSC70 by competing the co-chaperones.;induction:Constitutively synthesized.;PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.;similarity:Belongs to the heat shock protein 70 family.;subcellular location:Translocates rapidly from the cytoplasm to the nuclei, and especially to the nucleoli, upon heat shock. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.;subunit:Interacts with HSPH1/HSP105 . Interacts with IRAK1BP1 (By similarity). Interacts with PACRG and TSC2 . Interacts with SV40 VP1.;tissue specificity:Ubiquitous.;



Background


This gene encodes a member of the heat shock protein 70 family, which contains both heat-inducible and constitutively expressed members. This protein belongs to the latter group, which are also referred to as heat-shock cognate proteins. It functions as a chaperone, and binds to nascent polypeptides to facilitate correct folding. It also functions as an ATPase in the disassembly of clathrin-coated vesicles during transport of membrane components through the cell. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011],

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-HSC70 antibody. The HRP-conjugated Goat anti-Rabbit IgG (H + L) antibody was used to detect the antibody. Lane 1: HeLa Lane 2: Jurkat Lane 3: C6 Lane 4: NIH-3T3 Predicted band size: 71kDa Observed band size: 71kDa



Human breast carcinoma was stained with anti-HSC70 Rabbit antibody

