



Alkaline Phosphatase Rabbit mAb

Catalog No	YP-rAb-17179
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB,IHC,IF,IP,ELISA
Gene Name	ALPL
Protein Name	Alkaline phosphatase tissue-nonspecific isozyme
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	IHC 1:200-1:1000; WB 1:2000-1:10000; IF 1:200-1:1000; ELISA 1:5000-1:20000; IP 1:50-1:200; Note: For IHC, we suggest antigen retrieval with TE buffer pH 9.0
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	ALPL;Alkaline phosphatase tissue-nonspecific isozyme;ALPL ; Alkaline phosphatase ; tissue-nonspecific isozyme ; AP-TNAP ; TNSALP ; Alkaline phosphatase liver/bone/kidney isozyme
Observed Band	80kD
Calculated Molecular Weight	57kD
Cell Pathway	Cell membrane ; Lipid-anchor, GPI-anchor . Extracellular vesicle membrane ; Lipid-anchor, GPI-anchor . Mitochondrion membrane ; Lipid-anchor, GPI-anchor . Mitochondrion intermembrane space . Localizes to special class of extracellular vesicles, named matrix vesicles (MVs), which are released by osteogenic cells. Localizes to the mitochondria of thermogenic fat cells: tethered to mitochondria l membranes via a GPI-anchor and probably resides in the mitochondrion intermembrane space. .
Tissue Specificity	Brain,Cerebellum,Liver,Lymphoma,Osteosarcoma,Peripheral nerve,Semin
Function	Catalytic activity:A phosphate monoester + H(2)O = an alcohol + phosphate.,cofactor: Binds 1 magnesium ion.,cofactor: Binds 2 zinc ions.,Disease:Defects in ALPL are a cause of hypophosphatasia adult type (hypophosphatasia) [MIM:146300].,Disease:Defects in ALPL are a cause of hypophosphatasia childhood (hypophosphatasia) [MIM:241510].,Disease:Defects in ALPL are a cause of hypophosphatasia infantile (hypophosphatasia) [MIM:241500]; an inherited metabolic bone disease characterized by defective



skeletal mineralization. Four hypophosphatasia forms are distinguished, depending on the age of onset: perinatal, infantile, childhood and adult type. The perinatal form is the most severe and is almost always fatal. Patients with only premature loss of deciduous teeth, but with no bone disease are regarded as having odontohypophosphatasia (odonto).,Function:This isozyme may play a role in skeletal mineralization.,miscellaneous:In most mammals there are four different isozymes: placental, placental-like, intestinal and tissue non-specific (liver/bone/kidney).,online information:Alkaline phosphatase entry,online information:Tissue nonspecific alkaline phosphatase gene mutations database,PTM:Glycosylated.,similarity:Belongs to the alkaline phosphatase family.,subunit:Homodimer.,

Background


This gene encodes a member of the alkaline phosphatase family of proteins. There are at least four distinct but related alkaline phosphatases: intestinal, placental, placental-like, and liver/bone/kidney (tissue non-specific). The first three are located together on chromosome 2, while the tissue non-specific form is located on chromosome 1. The product of this gene is a membrane bound glycosylated enzyme that is not expressed in any particular tissue and is, therefore, referred to as the tissue-nonspecific form of the enzyme. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate the mature enzyme. This enzyme may play a role in bone mineralization. Mutations in this gene have been linked to hypophosphatasia, a disorder that is characterized by hypercalcemia and skeletal defects. [provided by RefSeq, Oct 2015]

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-Alkaline Phosphatase antibody. The HRP-conjugated Goat anti-Rabbit IgG (H + L) antibody was used to detect the antibody. Lane 1: HEK293 Lane 2: Mouse kidney Predicted band size: 57kDa Observed band size: 80kDa

