





DBX2 Rabbit mAb

Catalog No	YP-rAb-17420
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	DBX2
Protein Name	Homeobox protein DBX2
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-1:5000;
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	DBX2;Homeobox protein DBX2;Developing brain homeobox protein 2
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	similarity:Belongs to the H2.0 homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,
Background	similarity:Belongs to the H2.0 homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.





Western blot analysis of lysates from HepG2 cell, primary antibody was diluted at 1:1000, 4 ° over night, Dylight 800 secondary antibody was diluted at 1:10000, 37 ° 1 hour.



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4 ° over night

