




CECR5 Rabbit mAb

Catalog No	YP-rAb-17660
Isotype	IgG
Reactivity	Human
Applications	WB,IF,ICC,IP
Gene Name	HDHD5;CECR5
Protein Name	Haloacid dehalogenase-like hydrolase domain-containing 5;Cat eye syndrome critical region protein 5;
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-5000; ICC/IF 1:100-300; IP 1:50-100
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	HDHD5;Haloacid dehalogenase-like hydrolase domain-containing 5;Cat eye syndrome critical region protein 5;HDHD5 ; CECR5 ; Haloacid dehalogenase-like hydrolase domain-containing 5 ; Cat eye syndrome critical region protein 5 ;CECR5
Observed Band	47kD
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Widely expressed.
Function	
Background	Predicted to be involved in glycerophospholipid biosynthetic process. Located in mitochondrion. [provided by Alliance of Genome Resources, Feb 2025]
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Western Blot analysis of L02 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-CECR5 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

