




C3orf38 Rabbit mAb

Catalog No	YP-rAb-17662
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	C3orf38
Protein Name	Uncharacterized protein C3orf38
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-5000
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	C3orf38;Uncharacterized protein C3orf38;C3orf38 ; Uncharacterized protein C3orf38
Observed Band	36kD
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	May be involved in apoptosis regulation. {ECO:0000269 PubMed:17464193}.
Background	Involved in positive regulation of apoptotic process. Located in nucleus. [provided by Alliance of Genome Resources, Feb 2025]
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Western Blot analysis of L02 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-C3orf38 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

