




C2orf69 Rabbit mAb

Catalog No	YP-rAb-17663
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	C2orf69
Protein Name	Mitochondrial protein C2orf69
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-5000
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	C2orf69;Mitochondrial protein C2orf69;C2orf69 ; Mitochondrial protein C2orf69
Observed Band	42kD
Calculated Molecular Weight	
Cell Pathway	Mitochondrion matrix.
Tissue Specificity	
Function	May play a role in the respiratory chain. {ECO:0000269 PubMed:33945503}.
Background	Involved in oxidative phosphorylation. Located in mitochondrion. Implicated in combined oxidative phosphorylation deficiency 53. [provided by Alliance of Genome Resources, Feb 2025]
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.





Western Blot analysis of K-562 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-C2orf69 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

