




NUP62CL Rabbit mAb

Catalog No	YP-rAb-17664
Isotype	IgG
Reactivity	Human
Applications	WB,IF,ICC,FC,IP
Gene Name	NUP62CL;NUP62L
Protein Name	Nucleoporin-62 C-terminal-like protein
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-5000; ICC/IF 1:100-300; FC 1:100-300; IP 1:50-100
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	NUP62CL;Nucleoporin-62 C-terminal-like protein;NUP62CL ; NUP62L ; Nucleoporin-62 C-terminal-like protein ;NUP62L
Observed Band	20kD
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	This gene encodes a protein containing a domain found in nucleoporins which are glycoproteins found in nuclear pore complexes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2011]
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Western Blot analysis of HeLa whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-NUP62CL rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

