



RBM43 Rabbit mAb

Catalog No	YP-rAb-17666
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	RBM43;C2orf38
Protein Name	RNA-binding protein 43;RNA-binding motif protein 43;
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-5000
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	RBM43;RNA-binding protein 43;RNA-binding motif protein 43;RBM43 ; C2orf38 ; RNA-binding protein 43 ; RNA-binding motif protein 43 ;C2orf38
Observed Band	39kD
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	Predicted to enable RNA binding activity. [provided by Alliance of Genome Resources, Feb 2025]
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



A549

kDa

200-

140-

100-

70-

55-

40-


35-

25-

20-

15-

10-



RBM43

Western Blot analysis of A549 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-RBM43 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

