




BTNL3 Rabbit mAb

Catalog No	YP-rAb-17675
Isotype	IgG
Reactivity	Human
Applications	WB,FC
Gene Name	BTNL3;BTNL3;COLF4100;UNQ744/PRO1472
Protein Name	Butyrophilin-like protein 3;Butyrophilin-like receptor;
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-5000; FC 1:100-300
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	BTNL3;Butyrophilin-like protein 3;Butyrophilin-like receptor;BTNL3 ; BTNL3 ; COLF4100 ; UNQ744/PRO1472 ; Butyrophilin-like protein 3 ; Butyrophilin-like receptor ;BTNL3;COLF4100;UNQ744/PRO1472
Observed Band	51kD
Calculated Molecular Weight	
Cell Pathway	Membrane; Single-pass type I membrane protein.
Tissue Specificity	Expressed in small intestine, colon, testis, spleen, and leukocyte . {ECO:0000269 PubMed:10429365}.
Function	
Background	Predicted to enable signaling receptor binding activity. Predicted to be involved in T cell receptor signaling pathway and regulation of cytokine production. Predicted to be located in membrane. Predicted to be active in external side of plasma membrane. [provided by Alliance of Genome Resources, Feb 2025]
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Western Blot analysis of HeLa whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-BTNL3 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.





杭州臻优品生物科技有限公司

热销产品:

蛋白、一抗、抗体对、ELISA试剂盒、生化试剂盒
CCK8试剂盒、QPCR检测试剂盒

检测服务:

ELISA检测及定制服务 | 生化检测 | PCR、QPCR检测 | WB检测
ICO-IP检测 | 切片 | 染色 | 免疫组化 | 免疫荧光 | 透射电镜全套
| 宏基因组、转录组、基因组、蛋白组、代谢组测序



关注官网



关注客服