




TMEM192 Rabbit mAb

Catalog No	YP-rAb-17677
Isotype	IgG
Reactivity	Human
Applications	WB,FC
Gene Name	TMEM192
Protein Name	Transmembrane protein 192
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-5000; FC 1:100-300
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	TMEM192;Transmembrane protein 192;TMEM192 ; Transmembrane protein 192
Observed Band	30kD
Calculated Molecular Weight	
Cell Pathway	Lysosome membrane; Multi-pass membrane protein. Late endosome.
Tissue Specificity	Strongly expressed in kidney, liver, lung and pancreas . {ECO:0000269 PubMed:20370317}.
Function	
Background	Enables protein homodimerization activity. Located in several cellular components, including late endosome; lysosomal membrane; and perinuclear region of cytoplasm. [provided by Alliance of Genome Resources, Feb 2025]
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.





Western Blot analysis of HeLa whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-TMEM192 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

