




NBEAL2 Rabbit mAb

Catalog No	YP-rAb-17696
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	NBEAL2;KIAA0540;UNQ253/PRO290
Protein Name	Neurobeachin-like protein 2
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-5000
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	NBEAL2;Neurobeachin-like protein 2;NBEAL2 ; KIAA0540 ; UNQ253/PRO290 ; Neurobeachin-like protein 2 ;KIAA0540;UNQ253/PRO290
Observed Band	303kD
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum.
Tissue Specificity	Expressed in megakaryocytes. {ECO:0000269 PubMed:21765411}.
Function	Probably involved in thrombopoiesis. Plays a role in the development or secretion of alpha-granules, that contain several growth factors important for platelet biogenesis. {ECO:0000269 PubMed:21765411, ECO:0000269 PubMed:21765412}.
Background	The protein encoded by this gene contains a beige and Chediak-Higashi (BEACH) domain and multiple WD40 domains, and may play a role in megakaryocyte alpha-granule biogenesis. Mutations in this gene are a cause of gray platelet syndrome. [provided by RefSeq, Dec 2011]
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Western Blot analysis of HepG2 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-NBEAL2 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.





杭州臻优品生物科技有限公司

热销产品:

蛋白、一抗、抗体对、ELISA试剂盒、生化试剂盒
CCK8试剂盒、QPCR检测试剂盒

检测服务:

ELISA检测及定制服务 | 生化检测 | PCR、QPCR检测 | WB检测
| CO-IP检测 | 切片 | 染色 | 免疫组化 | 免疫荧光 | 透射电镜全套
| 宏基因组、转录组、基因组、蛋白组、代谢组测序



关注官网



关注客服