




TAGAP Rabbit mAb

Catalog No	YP-rAb-17698
Isotype	IgG
Reactivity	Human
Applications	WB,FC,IP
Gene Name	TAGAP;TAGAP1;FKSG15
Protein Name	T-cell activation Rho GTPase-activating protein;T-cell activation GTPase-activating protein;
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	WB 1:1000-5000; FC 1:100-300; IP 1:50-100
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	TAGAP;T-cell activation Rho GTPase-activating protein;T-cell activation GTPase-activating protein;TAGAP ; TAGAP1 ; FKSG15 ; T-cell activation Rho GTPase-activating protein ; T-cell activation GTPase-activating protein ;TAGAP1;FKSG15
Observed Band	80kD
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	May function as a GTPase-activating protein and may play important roles during T-cell activation. {ECO:0000269 PubMed:15177553}.
Background	This gene encodes a member of the Rho GTPase-activator protein superfamily. The encoded protein may function as a Rho GTPase-activating protein. Alterations in this gene may be associated with several diseases, including rheumatoid arthritis, celiac disease, and multiple sclerosis. Alternate splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2013],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Western Blot analysis of HEK-293 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-TAGAP rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

