



GAPDH Rabbit mAb

Catalog No	YP-rAb-17721
Isotype	IgG
Reactivity	Human,Mouse,Rat,Dog
Applications	WB,IHC,IF,ELISA
Gene Name	GAPDH
Protein Name	Glyceraldehyde-3-phosphate dehydrogenase
Purification Process	Protein A
Specificity	Endogenous
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Source	Monoclonal, Rabbit,IgG
Dilution	IHC 1:2000-1:5000; WB 1:10000-1:50000; IF 1:200-1:1000; ELISA 1:5000-1:20000; Note: For IHC, we suggest antigen retrieval with TE buffer pH 9.0
Concentration	0.5 mg/ml
Purity	≥90%
Storage Stability	-15 ° C to -25 ° C/1 year(Do not lower than -25 ° C)
Synonyms	GAPDH;Glyceraldehyde-3-phosphate dehydrogenase;GAPDH ; GAPD ; CDABP0047 ; OK/SW-cl.12 ; Glyceraldehyde-3-phosphate dehydrogenase ; GAPDH ; Peptidyl-cysteine S-nitrosylase GAPDH
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Cytoplasm, Nucleus
Tissue Specificity	Astrocytoma,Brain,Cajal-Retzius cell,Colon adenocarcinoma,Epitheliu
Function	Catalytic activity:D-glyceraldehyde 3-phosphate + phosphate + NAD(+) = 3-phospho-D-glyceroyl phosphate + NADH.,Function:Independent of its glycolytic activity it is also involved in membrane trafficking in the early secretory pathway.,online information:Glyceraldehyde 3-phosphate dehydrogenase entry,pathway:Carbohydrate degradation; glycolysis; pyruvate from D-glyceraldehyde 3-phosphate: step 1.,pathway:Carbohydrate degradation; glycolysis; pyruvate from D-glyceraldehyde 3-phosphate: step 1/5.,PTM:Reversible S-nitrosylation of Cys-152 inhibits enzymatic activity and increases endogenous ADP-ribosylation, which inhibits the enzyme in a non-reversible manner. The latter modification is more likely to be a pathophysiological event associated with inhibition of gluconeogenesis.,sequence



Caution: Differs quite extensively., similarity: Belongs to the glyceraldehyde-3-phosphate dehydrogenase family., subcellular location: Postnuclear and Perinuclear regions., subunit: Homotetramer., subunit: Homotetramer. Binds PRKCI.,

Background


glyceraldehyde-3-phosphate dehydrogenase (GAPDH) Homo sapiens This gene encodes a member of the glyceraldehyde-3-phosphate dehydrogenase protein family. The encoded protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions. The product of this gene catalyzes an important energy-yielding step in carbohydrate metabolism, the reversible oxidative phosphorylation of glyceraldehyde-3-phosphate in the presence of inorganic phosphate and nicotinamide adenine dinucleotide (NAD). The encoded protein has additionally been identified to have uracil DNA glycosylase activity in the nucleus. Also, this protein contains a peptide that has antimicrobial activity against *E. coli*, *P. aeruginosa*, and *C. albicans*. Studies of a similar protein in mouse have assigned a variety of additional functions including nitrosylation of nuclear proteins, the regulation of mRNA stability, and acting as a transferri

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Western Blot analysis using Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-GAPDH Rabbit mAb-YM8394 diluted at 1:20000. Secondary : Dylight 800, Goat Anti Rabbit IgG




Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-GAPDH antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: EL-4-B5 Lane 2: L6 Lane 3: JAR Lane 4: MDCK Lane 5: Cos-7 Predicted band size: 38kDa Observed band size: 38kDa






Human bladder carcinoma was stained with anti-GAPDH Rabbit antibody



Human liver was stained with anti-GAPDH Rabbit antibody



Loss of NDUFS1 promotes gastric cancer progression by activating the mitochondrial ROS-HIF1 α -FBLN5 signaling pathway. *BRITISH JOURNAL OF CANCER* Jin Zhou WB Human 1:5000 MKN45 cell, N87 cell, GES-1 cell, AGS cell, HGC-27 cell, KATO3 cell, SNU-1 cell

